

Issue date: 02-05-2025

Valid until: 02-05-2026

To whom it may concern

INFORMATION LETTER

Upon request of our authorized dealer on the territory of Ukraine, LLC Chemlaborreactive, we, YOUNG IN CHROMASS CO.,LTD, inform you that micro-TCD detector produced by our company is suitable for column flow up to 20 ml/min, can be used with packed, micropacked and capillary columns without any sensitivity or resolution loss.

We recommend to use micro-TCD detector type with packed columns for Natural Gas Analysis application according DSTU ISO 6974:2021.

02/05/2025

YOUNG IN CHROMASS CO.,LTD



JangWon Bae

YOUNGIN Chromass Co., Ltd.
60, Anyangcheondong-ro, Dongan-gu
Anyang-si, Gyeonggi-do, 14042, Korea



Дата видачі: 02-05-2025

Дійсний до: 02-05-2026

Тим, кого це стосується

ІНФОРМАЦІЙНИЙ ЛИСТ

На запит нашого авторизованого дилера на території України, ТОВ «Хімлаборреактив», повідомляємо, що мікро-ТCD-детектор виробництва нашої компанії підходить для потоку в колонці до 20 мл/хв, може використовуватися з набивними, мікронабивними та капілярними колонками без втрати чутливості чи роздільної здатності.

Ми рекомендуємо використовувати детектор типу micro-TCD з набивними колонками для аналізу природного газу відповідно до ДСТУ ISO 6974:2021.

02/05/2025

YOUNG IN CHROMASS CO.,LTD

Jang Won Bae

Jang Won Bae

YOUNGIN Chromass Co., Ltd.
60, Anyangcheondong-ro, Dongan-gu,
Anyang-si, Gyeonggi-do, 14042, Korea



ДИРЕКТОР
ЖОЛОБИЦЬКИЙ В.І.

